

SX

Autoclaves SX series
[SX-300] [SX-500] [SX-700]

Maximized chamber capacity at Minimum installation space



[SX-300]

[SX-500]

[SX-700]

NEW!

Easy operation using the "Top-open lid"



A top-open lid, which can be operated easily with a single hand or foot, is provided with the SX autoclaves. To open the lid, the lid lock is released by stepping on the foot pedal located at the bottom while pressing the lid down. The lid can be opened with minimal effort owing to mechanical assistance. In addition, the lid can be closed easily just by pressing it down.



Saving the installation space

The main body is very compact owing to the design which provides thorough installation space saving. In addition, installation space at the side of the main body which is required for models employing a slide-open lid, which is not necessary at SX autoclaves top-open lid. Even when several SX autoclaves are installed side-by-side, the necessary installation space is minimized.

Cooling fan provided as a standard feature

A rapid air-cooling function (vessel-cooling fan) permitting reduction of chamber cool down time is provided as a standard feature. This is most appropriate for lowering the temperature rapidly on completion of the cycle. The time required for lower the temperature is much reduced by employment of the cooling fan in comparison with natural cooling.

Sterilizing course selection

The optimal sterilizing course can be selected from among five courses, such as the liquid sterilizing course in addition to that for normal sterilization.

Liquid sterilizing course



This is most appropriate course setting for prevention of sudden sample boiling.

Sterilizing course



Course for normal sterilization.

Sterilizing-warming course



Prevents culture media from coagulating after sterilization.

Heating-warming course



Extremely convenient for dissolving culture media and warming.

Memory recall



Preferred operating conditions can be recalled up easily.



Easily viewable "work monitor"

The LED display shows working status for easy monitoring. The process being performed is indicated by an LED indication lamp blinking in red.



The illuminated display on the operating panel shown in the picture is different from the actual display at the time of operation.

Timer function

The operation starting time can be preset easily. By presetting the start of operation, otherwise wasted periods at night or early morning, can be used effectively.

OPTION

External sensor for the articles to be sterilized

Sensor for directly detecting the temperature of articles to be sterilized. More reliable sterilization monitor can be carried out.

Large indication lamp "Operated"

An easily viewable large indication "Operated" is provided at the upper section of the operating panel. Depending on the pressure status (normal pressure/pressure is applied), the displayed color changes.

Auto-variable exhaust speed

A function allowing the exhaust valve to open automatically after completion of sterilization is provided. The exhaust speed can be set to one of five levels (set to variable for liquid sterilization).

Lid interlock


A safety device (lid interlock mechanism) for locking the lid during operation is provided. A high level of safety is assured to prevent unanticipated accidents.

Data output

The temperature data, pressure data, etc. can be exported. The temperature can be recorded when the external recorder is connected.

Setting to 100 °C is possible.

Variable temperature setting such as to 100 °C as well as to 121 °C has been made possible. The temperature can be set within the range from 45 °C to 135 °C in steps of one degree (heating mode up to 104 °C).



Water level detector

A safety device (water level detector) for preventing dry heating is provided. The device is designed to detect the water level through the micro electrolytes of the sterilizing water.

External sensor and recorder

A temperature sensor and recorder independent of the autoclaves. The chamber temperature can be recorded.

Pressure fine adjustment function

When the balance between the temperature and pressure deviates during sterilization, fine exhausting is carried out automatically in order to adjust the best chamber condition.

Functions & Features

- Liquid sterilizing
- Heating
- Warming
- Check
- Last-run memory
- Built-in steam exhaust bottle
- Flat internal surface of the chamber
- Lid opening/closing detection mechanism
- Exhaust bottle detection mechanism
- Leakage breaker provided

Specifications

Model	SX-300	SX-500	SX-700
Operating temperature range	Sterilizing	105~135C (0.019MPa~0.212MPa)	
	Heating	45~104C (0~0.015MPa)	
	Warming	45~95C	
Maximum operating pressure	0.263MPa		0.25MPa
Temperature	Display	Digital	
Pressure gauge	Display	Analog	
	Display range	0~0.4MPa	
Heat source	1.5 kW electric heater	2.0 kW electric heater	3.0kW electric heater
Safety device	•Water level sensor •Current leakage breaker •Lid interlock •Over-heat prevention •Over-pressure prevention		
	•Open temperature sensor detection •Safety valve		
Time	Display	Digital	
	Display range	1 to 99 hours, 1 to 999 minutes (variable: 0:01 to 9:59/10 to 99)	
Chamber dimensions	Variable: 1 to 99 hours/ fixed to 4 hours		
	Effective internal volume: 36 L	Effective internal volume: 50 L	Effective internal volume: 69 L
	Internal volume of the chamber: 44 L	Internal volume of the chamber: 58 L	Internal volume of the chamber: 79 L
Chamber material	SUS304		
Dimensions (mm)	410W, 477D, 790H (With projection: 574 D)	410W, 477D, 970H (With projection: 574 D)	470W, 526D, 1003H (With projection: 625 D)
Weight	50kg	60kg	72kg
Rated voltage	120V AC		
Power input	13A/120V	17A/120V	—
	7A/230V	9A/230V	13A/230V
Required power supply	Single-phase 120 V AC (50/60 Hz) 15 A or more	Single-phase 120 V AC (50/60 Hz) 20A or more	—
	Single-phase 230 V AC (50/60 Hz) 10 A or more	Single-phase 230 V AC (50/60 Hz) 15A or more	Single-phase 230 V AC (50/60 Hz) 15A or more
Power consumption (caloric power)	1.5kW (1290kcal/h)	2kW (1720kcal/h)	3kW (2580kcal/h)
Environmental Conditions	When operating the autoclave, observe the environmental conditions given below.		
	Ambient temperature: 10 to 35C		
	Atmospheric pressure: 860 to 1060hPa		
	Relative humidity: 30 to 85%		
Accessories	Maximum gradient: 2'		
	Stainless basket (φ300 x 182 mm) 1	Stainless baskets (φ300 x 182 mm) 2	Stainless baskets (φ345 x 181mm) 2
Chamber bottom plate 1, Caster stoppers 4, Operation manual 1, Clear folder (for storing the operation manual) 1, Screw (for attaching the clear folder) 1, Warranty card 1, Customer card 1, Inspection sheet 1			

Accessories



Stainless Basket

The entire basket is made of punched metal and is designed for easy penetration of steam.



Stainless Bucket

Bucket without holes to prevent leakage of liquid from articles being sterilized.



Stainless Long Basket

Optimal for sterilizing the sterilizing bag. The structure of the lower section is without holes.



Stand for Testing Durham's Tubes

Convenient for sterilizing test tubes containing culture media. Can be used as a stand independently.



Sterilizing Bag

Sterilizing bags fit to the vessel dimensions of each model are available.

Sales Office:

TOMY DIGITAL BIOLOGY CO., LTD.
3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan
e-mail : info@digital-biology.co.jp
URL : http://www.digital-biology.co.jp
phone : +81-3-5971-8160 fax : +81-3-3970-6036

Manufacturer:

TOMY KOGYO CO., LTD.
3-14-17 Tagara, Nerima-ku, Tokyo 179-0073, Japan

TOMY TECH U.S.A., INC.

40479 Encyclopedia Circle, Fremont, California 94538, U.S.A.
e-mail : info@tomytech.com
URL : http://www.tomytech.com
phone : 510-440-1976 fax : 510-440-1975
Toll-Free US & Canada : 800-545-TOMY

TOMY SEIKO CO., LTD.

All TOMY products have a limited one-year warranty. Specifications are subject to change according to product advancement. Tomy and Digital Biology is registered trademark of Tomy Seiko Co., Ltd. and Tomy Digital Biology Co., Ltd. Copyright 2002, Tomy Seiko and its subsidiaries. Printed in Japan.

8K0120/3A

SX

Autoclaves SX series

[SX-300] [SX-500] [SX-700]



TOMY

SPEEDY AUTOCLAVES SX SERIES



 **Digital Biology®**